





























BRÁQUETE METÁLICO ADVANCED SERIES

PRESCRIÇÃO MBT 0,022"

SLOT 0,022" x 0,028" (0,56 x 0,71mm)

	Torque	-7°	-7°	0°	10°	17°	17°	10°	0°	-7°	-7°
	Ângulo	0°	0°	0°	8°	4°	4°	8°	0°	0°	0°
	Gancho	✓	✓	✓	-	-	-	-	✓	✓	✓
	Largura (mm)	2,8	2,8	2,8	2,8	3,1	3,1	2,8	2,8	2,8	2,8
		$\frac{5}{15}$	$\frac{4}{14}$	$\frac{3}{13}$	$\frac{2}{12}$	$\frac{1}{11}$	$\frac{1}{21}$	$\frac{2}{22}$	$\frac{3}{23}$	$\frac{4}{24}$	$\frac{5}{25}$
											
	Maxilar	Superior Direito Inferior Direito					Superior Esquerdo Inferior Esquerdo				
	Mandibular										
		$\frac{45}{5}$	$\frac{44}{4}$	$\frac{43}{3}$	$\frac{42}{2}$	$\frac{41}{1}$	$\frac{31}{1}$	$\frac{32}{2}$	$\frac{33}{3}$	$\frac{34}{4}$	$\frac{35}{5}$
	Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
	Ângulo	2°	2°	3°	0°	0°	0°	0°	3°	2°	2°
	Gancho	✓	✓	✓	-	-	-	-	✓	✓	✓
	Largura (mm)	2,8	2,8	2,8	2,2	2,2	2,2	2,2	2,8	2,8	2,8